

Sensitive Analysis of Neutral N-Glycan using Anion Doped Liquid Matrix by Negative-Ion Matrix-Assisted Laser Desorption/Ionization Mass Spectrometry

O Takashi Nishikaze, Yuko Fukuyama, Shin-ichirou Kawabata, Koichi Tanaka

Koichi Tanaka Laboratory of Advanced Science and Technology, Shimadzu Corporation, Kyoto, Japan

1. Overview

- ✓ Combinations of five liquid matrices and nine anions were evaluated for sensitive analyses of underivatized neutral *N*-glycans in negative-ion MALDI-MS.
- ✓ Among tested, G_3CA exhibited outstanding performances for all the anions.
- $\checkmark\,$ Low-femtomole or sub-femtomole detection limits could be achieved using NO_3^ or BF_4^ doped G_3CA.
- $\checkmark\,$ A combination of G_3CA and NO_3^- is the most suitable for sensitive analyses, especially for MS^2 measurement.
- ✓ A combination of G₃CA and BF₄⁻ is a possible alternative for superior sensitivity in MS¹ detection and medium sensitivity in MS² measurement.



2. Introduction

Negative-ion fragmentation of *N*-glycans has been proven to be more informative than that of positive-ion. It defines structural features such as the specific composition of the two antennae, the location of fucose, and the presence or absence of bisecting GlcNAc. In general, negative ionization of neutral *N*-glycans is much more difficult than positive ionization because of lack of acidic groups. The most promising approach is to ionize them as anionic adducts. However, anionic adduct formation of neutral *N*-glycans by MALDI-MS remains a challenging task, and the detection limit in negative-ion mode is merely at the sub-picomole level. Here, we demonstrate that G₃CA acts as an effective matrix for sensitive MALDI analysis of neutral *N*-Glycan in negative-ion mode.



3.3. Preparation for liquid matrices

Matrix (DHB, THAP, HABA, or CA) + counter cation (BA or TMG) in MeOH ↓ SpeedVac (overnight)

↓ Dissolve in 50% ACN at a concentration of 100mg/mL 20-50 : 1 Dilute fivefold with 50% ACN containing 2.5mM anions (mol/mol)

3.4. Sample preparation and mass spectrometry

- ✓ On target mix (Sample solution 0.5µL + matrix solution 0.5µL)
- ✓ µFocus MALDI plate[™] 700µm (Hudson Surface Technology)
- ✓ MALDI-QIT-TOF-MS (AXIMA Resonance,[™] Shimadzu/Kratos)

4. Results and Discussion

4.1 Evaluation of anion-doped liquid matrices

Combinations of five liquid matrices and nine anions were used to ionize *N*-glycans as anionic adducts, and their performances for sensitive analyses were evaluated. Representative spectra were shown in Figure 1 (below).

			The results obtained in this study could be interpreted		
Anion	Anion Nominal GB mass (kJ/mol)		based on the gas-phase basicity (GB) of anions.		
PF_6^-	145	1157		✓ Strong signals of anion-adducted NA2 glycans in combinations with any liquid	
BF_4^-	87	1204	Anions		
HSO4	97	1265	Low GB	matrices	used.
ľ	127	1294		✓ No or ne	egligible fragment ions
NO3	62	1330		✓ DHBB, GTHAP, and G ₂ CHCALittle or no signals of adducts	
Br	79, 81	1332	Anions with High GB		
SCN	58	1343		✓ G ₂ HABA	Weak signals of adducts
H ₂ PO ₄	97	1351		and exte	ensive fragment ions
Cľ	35, 37	1373 🦯)	✓ G ₃ CA	Strong signals of adducts
RA(%) DHE 50 GB 1322 0 GTH 50 G2C 0 G2H	BF ₄ - do BB of [DHB-H]-; kJ/mol) AP of [THAP-H]- kJ/mol) HCA	1727.6 1727.6 1727.6 1727.6	NO ₃ -	doped	GB was important parameter. Because ✓ Anionic adducts of glycans can be considered as proton-bound mixed dimers of anions, [M-H]····H [*] ···[A] [*] . ✓ Anions with high GB are easily neutralized by matrix molecules. Matrix + [A] [*]
G ₃ C	A	1727.6		1702.6 1702.6 1600 1800 200 m/2	→ [Matrix-H] + HA Figure 1. Negative-ion mass spectra of NA2 glycan at 100fmol using various liquid matrices and anions as dopants.





Since fragmentation of anionic adducts was initiated by proton abstraction from *N*-glycan to anions, the GB of anions is a crucial parameter for successful MS² analysis. Anions with low GB are unable to abstract hydroxylic protons from *N*-glycans, resulting in poor or no glycan fragment ions in the spectra.



100fmol NA2 + low GB anions

 $(PF_6^-, HSO_4^-, and I^-)$

✓ Poor MS² spectra

✓ Need extensive spectral accumulation

✓ Exception; BF₄- (GB 1204 kJ/mol)

 ✓ Glycan fragments with moderate efficiency
C_{3a or β}^{-3A}/₋₄



100fmol NA2 + high GB anions

(NO3-, Br-, SCN-, H2PO4-, and Cl-)

 \checkmark MS² spectra with high quality

 $\checkmark \rm NO_3^-$ was the best

MS

- \checkmark Strong signals of precursor ions
- ✓ Efficient fragmentation





Figure 4. Negative-ion MS² spectra of the anion-adducted NA2 glycans at 100fmol.

The I- adduct did not exhibit any informative signals in the MS² spectrum, whereas the BF. adduct produced MS² spectra with an acceptable s/n ratio MS² of the BF₄⁻ adduct formed product ions of glycans more efficiently than did the I- adduct. Using NO3-, the MS² spectrum was successfully obtained with a sufficient s/n ratio Since the NO2- adduct is less sensitive than the I or BF4 adduct in MS¹ analyses (Figure 2.), these results clearly indicate the effectiveness of NO3⁻ adduct for sensitive MS² measurement.

5. Conclusions

✓ A combination of G₃CA and NO₃⁻ is the most suitable for sensitive analyses, especially for MS² measurement of neutral *N*-glycans in negative-ion MALDI-MS.

✓ A combination of G_3CA and BF_4^- is a possible alternative for superior sensitivity in MS¹ detection and medium sensitivity in MS² analysis.

6. Acknowledgments

This research is supported by the Japan Society for the Promotion of Science (JSPS) through its Funding Program for World-Leading Innovative R&D on Science and Technology (FIRST Program).