


1. Overview.

- Alkylated trihydroxyacetophenone (ATHAP) as a novel matrix for hydrophobic peptides was reported.
- ATHAP is a 2,4,6-trihydroxyacetophenone derivative incorporating a hydrophobic alkyl chain on the acyl group which was expected to have affinity for hydrophobic peptides.
- As a result, ATHAP increased the sensitivity of hydrophobic peptides by 10-fold and decreased the sensitivity of hydrophilic peptides compared with α -cyano-4-hydroxycinnamic acid (CHCA).
- The peptides were detected throughout the entire matrix/analyte dried spot using ATHAP.

2. Introduction.

- Hydrophobic peptides are difficult to be detected in MALDI-MS because the detection hindrance appertains to hydrophilic peptides of conventional matrices.
- Recently, we reported alkylated dihydroxybenzoic acid (ADHB) as a matrix additive for hydrophobic peptides. However, it still remains the following issues:
 - Hydrophobic peptide ions were detected in the rim of matrix/analyte dried spot, which has difficulty in finding the "sweet spot".
 - ADHB was an additive, thus unavailable without conventional matrices.
 - Hydrophilic peptide ions were also detected, which may limit the detection of hydrophobic peptide ions.
- To solve the issues, we launched a study to develop a novel matrix for hydrophobic peptides.



3. Methods.

3-1. Alkylated trihydroxyacetophenone (ATHAP).

- ATHAP incorporating a C8 acyl chain (Figure 1) was synthesized.
- ATHAP solution was prepared in 75% acetonitrile (ACN)/0.1% aqueous trifluoroacetic acid (TFA) (v/v) at 5 mg/mL.

3-2. Matrix solution.

- α -Cyano-4-hydroxycinnamic acid (CHCA) was purchased from LaserBio Labs.
- Trihydroxyacetophenone (THAP) was purchased from Sigma-Aldrich.
- Each matrix was dissolved in 50% ACN/0.1% aqueous TFA (v/v) at 10 mg/mL.

3-3. Analyte solution.


- The peptides were dissolved in 50% ACN/0.1% aqueous TFA at appropriate concentrations.

3-4. Sample preparation.

- The analyte solution (0.5 μ L) and the matrix solution (0.5 μ L) were mixed on a stainless-steel plate to be analyzed by MALDI-TOFMS (Scheme 1).

3-5. MALDI-MS.

- MALDI-TOFMS measurement was performed using AXIMA Performance (Shimadzu/Kratos, UK) mass spectrometer in linear, positive ion mode.



4. Results.

4-1. Sensitivity Increase for Hydrophobic Peptides by ATHAP.


Table 1. Detection limit of hydrophobic peptide.

| | humanin (fmol/well) |
|-----------|---------------------|
| C12-ATHAP | 10 |
| C10-ATHAP | 1 |
| C8-ATHAP | 1 |
| C6-ATHAP | 10 |
| THAP | 10 |
| CHCA | 10 |


The sensitivity for hydrophobic peptide was increased by 10-fold using C8- or C10-ATHAP compared with THAP or CHCA.

C8 or C10

4-2. Preferential Ionization of Hydrophobic Peptides by ATHAP.



4-3. Enhanced Detection of [M+H]⁺ throughout the Entire Spot by ATHAP.



4-4. Sensitivity Improvement by ATHAP and SSRCalc Hydrophobicity.

- The sensitivity improvement rate by ATHAP to CHCA was increased for hydrophobic peptides with higher SSRCalc Hydrophobicity and decreased for hydrophilic peptides with lower SSRCalc Hydrophobicity.

Table 2. Sensitivity improvement by ATHAP for peptides with SSRCalc Hydrophobicity from 5.2 to 54.8.

| no. | name | analyte | | |
|-----|------------------------------|------------------------|------------|--|
| | | SSRCalc Hydrophobicity | m/z (Ave.) | sensitivity improvement rate by ATHAP to CHCA (fold) |
| 1 | NF- κ B inhibitor | 54.8 | 2782.6 | 10 |
| 2 | OVA-BIP hybrid peptide | 50.2 | 2291.6 | 10 |
| 3 | humanin | 50.0 | 2688.3 | 10 |
| 4 | β -amyloid 22-42 | 42.4 | 2000.4 | 10 |
| 5 | catestatin | 38.1 | 2327.7 | 1 |
| 6 | ACTH 18-39 | 37.9 | 2466.7 | 1 |
| 7 | nocistatin | 29.8 | 1928.1 | 1 |
| 8 | neuropeptide S | 22.3 | 2188.5 | 0.1 |
| 9 | β -amyloid 1-16 | 18.2 | 1956.1 | 0.1 |
| 10 | β -amyloid 1-11 | 13.5 | 1326.3 | 0.001 |
| 11 | β -conglycinin 165-178 | 5.2 | 1848.8 | 0.0001 |

Detection limit using CHCA

Detection limit using ATHAP

Higher sensitivity using ATHAP than CHCA

Lower sensitivity using ATHAP than CHCA

References.

[1] Fukuyama, Y.; Tanimura, R.; Maeda, K.; Watanabe, M.; Kawabata, S.; Iwamoto, S.; Izumi, S.; Tanaka, K. *Anal. Chem.* **2012**, *84*, 4237-4243.

Acknowledgments.

This research is granted by the Japan Society for the Promotion of Science (JSPS) through the "Funding Program for World-Leading Innovative R&D on Science and Technology (FIRST Program)," initiated by the Council for Science and Technology Policy (CSTP).

4-5. Phosphorylase b Lys-C Digestion Fragments Analysis using ATHAP.




Table 3. The ion detection for phosphorylase b Lys-C digestion fragments and the SSRCalc Hydrophobicity.^a

| No | SSRCalc Hydrophobicity | m/z (Ave.) | ATHAP | CHCA |
|----|------------------------|------------|-------|------|
| 1 | 55.9 | 3602.2 | ** | - |
| 2 | 53.9 | 3823.5 | ** | - |
| 3 | 53.1 | 3890.3 | ** | - |
| 4 | 51.0 | 3823.5 | ** | - |
| 5 | 50.7 | 2198.6 | ** | + |
| 6 | 45.8 | 2155.6 | ** | + |
| 7 | 45.1 | 2742.0 | + | + |
| 8 | 42.8 | 2969.5 | + | + |
| 9 | 38.9 | 3504.9 | + | + |
| 10 | 35.1 | 1855.1 | ** | ** |
| 11 | 33.9 | 1657.0 | + | + |
| 12 | 33.6 | 2130.5 | ** | + |
| 13 | 33.4 | 2620.0 | ** | ** |
| 14 | 31.4 | 1610.9 | ** | ** |
| 15 | 31.4 | 1814.1 | ** | ** |
| 16 | 31.2 | 1526.8 | ** | ** |
| 17 | 30.9 | 2043.3 | ** | ** |
| 18 | 30.4 | 1942.3 | ** | ** |
| 19 | 29.9 | 2449.7 | - | + |
| 20 | 28.4 | 1304.7 | ** | ** |
| 21 | 25.8 | 1263.4 | - | ** |
| 22 | 24.7 | 1178.3 | - | ** |
| 23 | 24.5 | 1290.5 | + | ** |
| 24 | 18.9 | 1254.5 | + | ** |
| 25 | 9.3 | 1102.2 | - | ** |

Higher intensity using ATHAP than CHCA

Lower intensity using ATHAP than CHCA

Hydrophobic peptide ions, which could not be detected by CHCA, were detected by ATHAP.

- The sequence coverage was increased by using the both results of ATHAP and CHCA.
 - 30% with CHCA
 - 45% (or 49%) with ATHAP
 - 51% (or 55%) with both the results of ATHAP and CHCA

Higher intensity using ATHAP than CHCA.

Lower intensity using ATHAP than CHCA.

^a "*" indicates that the ions were detected with S/N ≥ 5 . "+" indicates that the ions were detected with S/N = 2-5, and "-" indicates that the ions were not detected.

5. Conclusions.

- A novel matrix ATHAP increased sensitivity of hydrophobic peptides by 10-fold.
- All the issues of ADHB were solved by ATHAP as follows:
 - The peptide ions were detected throughout the entire matrix/analyte dried spot.
 - ATHAP works as a matrix itself.
 - The detection of hydrophilic peptide ions was suppressed.