G SHIMADZU Improvement of MALDI-TOFMS Ion Source for Higher Resolving Power over Wider Mass Range

Kei Kodera¹, Toshinori Kobayashi¹, Sadanori Sekiya¹, Makoto Hazama¹, Shinichi Iwamoto¹, Koichi Tanaka¹ 1 Shimadzu Corporation. 1, Nishinokyo-Kuwabaracho Nakagyo-ku, Kyoto 604-8511, Japan

1. Introduction

Delaved pulse extraction is commonly used as a method which compensates for the initial kinetic energy scattering to obtain higher resolving power of time of flight mass spectrometer (TOEMS) with matrix-assisted laser desorption/ionization (MALDI). However, the effective mass range of this method is limited and not wide enough. We have previously developed a new method which expands the effective range by dividing ions according to their masses before applying pulse voltage. In this study, we improved the method to further expand mass range by non-linear field

2. Methods

2-1. Pre-extraction method



Figure 1 Extraction method

Figure 1 shows the schemes of a) conventional and that of b) pre-extraction delayed pulse method. The initial ion velocity is independent of the mass, which is a feature of MALDI. Therefore the spatial distribution during the delay time is also independent of the mass without the pre-extraction as in figure 1 a), and each mass needs to be supplied optimum energy by pulse field. That is, effective mass range by specific pulse is narrow. On the other hand, by the separation of preacceleration field, the effective mass range of pre-extraction method is wider (Figure 1 b)

2-2. Optimum potential distribution

Before the pulse voltage is applied, the electric field E₀ between sample plate and 1st electrode (extraction electrode) is constant. The potential $V_0(x)$ is expressed as $V_0(x) = V_c - E_0 x$ (1)

where V_s is the initial sample voltage and x is the distance from sample plate surface. The flight time t = 0 is defined as the time when ions are generated on the sample plate surface. When the initial velocity is v_0 , it may have a distribution from $v_0 - \Delta v_0$ to $v_0 + \Delta v_0$. Thus, the difference in kinetic energy ΔK between ions which have initial velocities $v_0 + \Delta v_0$ and $v_0 - \Delta v_0$ is



 $\Delta K = \frac{1}{2}m(v_0 + \Delta v_0)^2 - \frac{1}{2}m(v_0 - \Delta v_0)^2$

 $= 2mv_0\Delta v$

where m is the ion mass.

$$\Delta U = -2\Delta x \times q \frac{d}{dx} V_{\rho}(x) = -2q\Delta v_{0}t_{0}\frac{d}{dx}V_{\rho}(x).$$

When ΔU compensates for ΔK , the equation would be
 $\frac{d}{dx}V_{\rho}(x) = -\frac{m}{q}\frac{v_{0}}{t_{0}}$

from Eqs. (2) and (6), where
$$\Delta U = \Delta K$$
. And substituting Eq. (5) into Eq. (7) gives

$$\frac{d}{d} V_{p}(\mathbf{x}) = -\frac{E_{0} v d_{0}}{2} \frac{1}{1 - \frac{1}{2}}.$$
(8)

$$V_p(x) = -\frac{E_0 v_0 t_0}{2} \ln(x - v_0 t_0) + V_c.$$

When the potential at
$$t = t_0$$
 defined $V_1(x)$, we can write it as
 $V_1(x) = V_0(x) + V_p(x) = V_s - E_0 \left\{ x + \frac{v_0 t_0}{2} \ln(x - v_0 t_0) \right\}$

from Eq. (9), where $V_s' = V_s + V_c$.

So at

Initial

line

Eq. (10) does not include mass-to-charge ratio m/q. Consequently, by applying pulse voltage which satisfies Eq. (10), we can compensate any mass ion for initial kinetic energy difference. However it doesn't work for time focusing since ions fly at same velocity in the drift region. Though further compensation is required for the ideal time focusing, it is clear that the compensation is based on V1(x) and the form



Figure 2 Ideal and approximate potential





(2)

(6)

(7)

(9)

(10)







Table 1 Parameters of each electrodes in simulations

Figure 3 Schematic representation

of linear TOF setun

	Sample	Sample plate			action de		2 nd extraction	acceleration
	Base (kV)	Pulse (V)	Delay (ns)	Base (kV)	Pulse (V)	Delay (ns)	(kV)	electrode (V)
Convention: method	al 18.0	1200	600	18.66	-	-	5.0	0
Method 1	18.0	800	450	17.2	-		5.0	0
Method 2	18.0	1540	730	17.2	830	730	16.0	0



4. Conclusions





ASMS2013 ThP06-088



3-2. Experimental Results

separated when method 2 was used.

Table 2 Parameters of each electrodes in experiments Delay

(ns)

730 18.66

400

1st extraction

Base Pulse

electrode

(kV) (V) (ns)

17.5

2nd extraction

ectrode

(kV)

4.7

6.0

Delay

acceleration

ectrode

(V)

0

0









Figure 8 Resolving power of PMMA peaks with each method

Acknowledgements

This research is granted by the Japan Society for the Promotion of Science (JSPS) through the "Funding Program for World-Leading Innovative R&D on Science and Technology (FIRST Program)," initiated by the Council for Science and Technology Policy (CSTP).