

# シアル酸へのin-gelメチルアミド化を用いた血清IgGのN型糖鎖プロファイリング N-glycan profiling of serum IgG by in-gel methylation for sialic acids

○金城 薫・濱名 周子・西風 隆司・川畑 慎一郎・田中 耕一 株式会社島津製作所 田中最先端研究所

○Kaoru Kaneshiro, Chikako Hamana, Takashi Nishikaze, Shin-ichirou Kawabata, Koichi Tanaka Koichi Tanaka Laboratory of Advanced Science and Technology, Shimadzu Corporation

## 1. Introduction

Matrix-assisted laser desorption/ionization (MALDI) mass spectrometry (MS) is an useful tool for the structural analysis of glycoproteins. MS1 and MS2 analyses of *N*-glycans released from glycoproteins are performed to obtain both glycan profiles and detailed structure of glycans. We previously reported the on-target AQ-labeling technique using a liquid matrix as a highly sensitive method for *N*-glycan profiling.<sup>1,2</sup> Recently we also reported that AQ-labeled *N*-glycans exhibit simple and informative negative-ion CID spectra with product ions due to cross-ring cleavages of the chitobiose core and ions specific to two antennae (D and E ions).<sup>3</sup>

The major limitation of current *N*-glycan profiling methods is the partial loss of sialic acid residues in MS<sup>1</sup> experiment. So far, to suppress the dissociation, many derivatization methods for sialic acid have been proposed such as methylesterification, methylation, and amidation. However, most of these methods require complicated purification procedures prior to MS.

Here, we developed an in-gel derivatization method for sialylated glycoprotein following SDS-PAGE, a common method of electrophoresis applied to isolate a purposive protein from complex biological samples. Among tested, the combination of the novel in-gel methylation method and the AQ-labeling method was determined to be best for the structural analysis of sialylated *N*-glycans.

- (1) Kaneshiro, K. et al. *Anal. Chem.* **2011**, *10*, 3663-7.
- (2) Kaneshiro, K. et al. *Anal. Chem.* **2012**, *16*, 7146-51.
- (3) Nishikaze, T. et al. *Anal. Chem.* **2012**, *21*, 9453-61.


## 2. Experimental Section

### 2-1. Materials

- ✓ Transferrin Human / Sigma-Aldrich
- ✓ Serum – Normal human – Male (single donor) / Tokyo Future Style, Inc.
- ✓ Serum – Normal human – Female (single donor) / Tokyo Future Style, Inc.
- ✓ Protein G Sepharose™4 Fast Flow / GW Healthcare Life Sciences
- ✓ NuPAGE® LDS Sample Buffer (x4) / Invitrogen
- ✓ NuPAGE® Novex 10% Bis-Tris Gel / Invitrogen
- ✓ Methylamine hydrochloride / Sigma-Aldrich
- ✓ 4-(4,6-Dimethoxy-1,3,5-triazin-2-yl)-4-methylmorpholinium chloride n-Hydrate (DMT-MM) / Wako

### 2-2. Isolation of IgG from serum


#### 【 Schematic diagram of IgG Isolation from serum 】




### 2-3. MALDI-MS analysis of sialoglycans by in-gel methylation

#### followed by AQ-labeling method

#### 【 Schematic diagram of in-gel methylation and AQ-labeling of sialoglycans 】



#### 【 3-AQ/CA liquid matrix<sup>4</sup> 】



(4) Fukuyama, Y et al. *59<sup>th</sup> ASMS Conference on Mass Spectrometry and allied topics*, June 5-9, 2011, ThP583

## 3. Results and Discussion

### 3-1. In-gel derivatization efficiency of sialic acids on glycoprotein

Conventional in-solution derivatization methods for sialoglycan were applied to in-gel derivatization for sialylated glycoprotein.

- ◆ The choice of solvent was crucial. In-gel derivatization reaction could only be performed in water (not in MeOH or DMSO).





Figure 1 Negative-ion MS spectra of in-gel derivatized transferrin *N*-glycan: (a) methylesterification (Me); (b) derivatization with acetohydrazide (AH); (c) amidation (Am); (d), (e) methylamidation (MA).

### 3-3. Optimum condition of in-gel methylation

In-gel methylation condition was tested at various times and temperatures to determine the optimum protocol.

- ◆ As the heating temperature rises, the reaction of methylation is promoted.
- ◆ The optimal results were obtained at 80 °C for 18 h.



### 3-5. Analysis of structure of in-gel methylated IgG glycan

The detail of methylated IgG glycan structure were analyzed by negative ion CID-MS.

- ◆ D ion was clearly detected in MS<sup>2</sup> spectra.
- ◆ MS<sup>2</sup> were successfully obtained from possible low intensity ions from MS<sup>1</sup> profiles.

### 3-2. Suitability for negative-ion CID-MS

- ◆ Derivatization with acetohydrazide produced multiple cleavages and faint D ion in MS<sup>2</sup> spectra.
- ◆ Methyl ester form of sialoglycan showed extensive loss of MeOH molecules, interfering with correct MS<sup>2</sup> spectral interpretation.
- ◆ Methylamidated sialoglycans exhibited simple and clear MS<sup>2</sup> spectra.




Figure 2 Negative-ion MS<sup>2</sup> spectra of in-gel derivatized A2 glycans: (a) derivatization with acetohydrazide; (b) methyl esterification; (c) methylamidation.

### 3-4. N-glycan profiles of IgG isolated from serum

Isolation of IgG from serum ⇒ in-gel methylation ⇒ *N*-glycan profiling

- ◆ IgG sialoglycans were completely methylated.
- ◆ We successfully profiled *N*-glycans of IgG isolated from 1 µL of human serum (male and female).




Figure 4 Negative-ion MS spectra of methylated *N*-glycan profiles of IgG isolated from human serum: (a) male; (b) female.

Table 1 Contents of each glycan relative to total 14 glycans ion intensities.

Glycan No.	contents (%)	
	Male	Female
G1	1.32	1.50
G2	31.63	29.50
G3	2.54	1.71
G4	39.10	40.36
G5	15.71	17.02
G6	2.43	3.04
G7	1.19	1.26
G8	0.81	0.80
G9	0.87	0.71
G10	3.13	2.98
G11	0.24	0.24
G12	0.58	0.49
G13	0.29	0.25
G14	0.15	0.13
Total	100.00	100.00




Figure 5 Negative-ion MS<sup>2</sup> spectra of in-gel methylated IgG *N*-glycans isolated from male serum. The selected precursor ions were (a) *m/z* 1847.6 (G4) and (b) *m/z* 2821.1 (G14)

## Conclusions

- ✓ In-gel complete derivatization for sialic acids of IgG was accomplished with an optimal condition of methylation using methylamine hydrochloride and DMT-MM (condensing agent).
- ✓ This novel method requires only washing gel pieces with water, while the conventional in-solution derivatization methods for sialylated glycans involve complicated purification processes to remove reagents and salts. ⇒ Preventing loss of sample caused by time-consuming purification processes (e.g. solid-phase extraction).
- ✓ *N*-glycan profiling was obtained without a loss of sialic acid residues in MS<sup>1</sup> analysis.
- ✓ In-gel methylation followed by on-target AQ-labeling method was highly suitable for the structural characterization of IgG *N*-glycans.

## Acknowledgement

This research is granted by the Japan Society for the Promotion of Science (JSPS) through the "Funding Program for World-Leading Innovative R&D on Science and Technology (FIRST Program)," initiated by the Council for Science and Technology Policy (CSTP).