

# 疎水性ペプチドに適したMALDIマトリックスalkylated trihydroxyacetophenone Alkylated Trihydroxyacetophenone as a MALDI matrix for Hydrophobic Peptides

○ 福山裕子<sup>1</sup>・中島ちひろ<sup>1</sup>・古市圭子<sup>2</sup>・谷口謙一<sup>1</sup>・川畑慎一郎<sup>1</sup>・泉俊輔<sup>2</sup>・田中耕一<sup>1</sup> <sup>1</sup>島津製作所 田中最先端研究所・<sup>2</sup> 広島大学 理学研究科  
○ Yuko Fukuyama<sup>1</sup>, Chihiro Nakajima<sup>1</sup>, Keiko Furuichi<sup>2</sup>, Kenichi Taniguchi<sup>1</sup>, Shin-ichirou Kawabata<sup>1</sup>, Shunsuke Izumi<sup>2</sup>, Koichi Tanaka<sup>1</sup>  
<sup>1</sup>Koichi Tanaka Laboratory of Advanced Science and Technology, Shimadzu Corporation, <sup>2</sup>Department of Mathematical and Life Sciences, Hiroshima University




SHIMADZU

ms<sup>3</sup>d  
FIRST Program

## 1. Introduction


- Hydrophobic peptides are difficult to be detected in MALDI-MS because the detection hindrance appertains to hydrophilic peptides of conventional matrices.
- Recently, we reported alkylated dihydroxybenzoic acid (ADHB) as a matrix additive for hydrophobic peptides. But, it still remains the following issues:

- Hydrophobic peptide ions were detected in the rim of matrix/analyte dried spot, which has **difficulty in finding the "sweet spot"**.
- ADHB was an additive, thus **unavailable without conventional matrices**.
- Hydrophilic peptide ions were also detected**, which may limit the detection of hydrophobic peptide ions. **A matrix additive ADHB<sup>[1]</sup>**

A matrix additive ADHB<sup>[1]</sup>

- To solve the issues, we launched a study to develop a novel matrix for hydrophobic peptides.<sup>[2]</sup>

## 2. Experimental Section

Fig 1. A novel matrix **ATHAP**

### 2-1. Alkylated trihydroxyacetophenone (ATHAP).

- ATHAP incorporating a C8 (or C6, C10, C12) acyl chain (Figure 1) was synthesized.
- ATHAP solution was prepared in 75% acetonitrile (ACN)/0.1% aqueous trifluoroacetic acid (TFA) (v/v) at 5 mg/mL.

### 2-2. Matrix solution.

- $\alpha$ -Cyano-4-hydroxycinnamic acid (CHCA) was purchased from LaserBio Labs.
- Trihydroxyacetophenone (THAP) was purchased from Sigma-Aldrich.
- Each matrix was dissolved in 50% ACN/0.1% aqueous TFA (v/v) at 10 mg/mL.

### 2-3. Analyte solution.


- The peptides were dissolved in 50% ACN/0.1% aqueous TFA (v/v) at appropriate concentrations.

### 2-4. Sample preparation.

- The analyte solution (0.5  $\mu$ L) and the matrix solution (0.5  $\mu$ L) were mixed on a stainless-steel plate to be analyzed by MALDI-TOFMS (Scheme 1).

### 2-5. MALDI-MS.

- MALDI-TOFMS measurement was performed using AXIMA Performance (Shimadzu/Kratos, UK) mass spectrometer in linear, positive ion mode.



Scheme 1. Sample preparation

## 3. Results and Discussion

### 3-1. ATHAP for Hydrophobic peptides




Fig 2. Mass spectra for 100 fmol equimolar mixture of hydrophobic peptide (humanin) and hydrophilic peptide ( $\beta$ -amyloid 1-11). ND indicates that the ions were not detected.




Fig 3. MS imaging of  $[M+H]^+$  for 100 fmol humanin.

- Hydrophobic peptide was preferentially ionized using ATHAP.
- The detection of hydrophilic peptide ions was suppressed using ATHAP.
- Hydrophobic peptide ions were detected throughout the entire matrix/analyte dried spot using ATHAP.  $\rightarrow$  Raster scanning was used in all the experiments.

Table 1. Sensitivity improvement by ATHAP for peptides

no.	name	analytes		sensitivity improvement rate by ATHAP to CHCA (fold)	Detection limit using CHCA = Detection limit using ATHAP
		SSRCalc Hydrophobicity <sup>[3]</sup>	m/z (Ave.)		
1	NF- $\kappa$ B inhibitor	54.8	2782.6	10	Higher sensitivity with ATHAP
2	OVA-BIP hybrid peptide	50.2	2291.6	10	
3	humanin	50.0	2688.3	10	
4	$\beta$ -amyloid 22-42	42.4	2000.4	10	
5	catestatin	38.1	2327.7	1	Lower sensitivity with ATHAP
6	ACTH 18-39	37.9	2466.7	1	
7	nocistatin	29.8	1928.1	1	
8	neuropeptide S	22.3	2188.5	0.1	
9	$\beta$ -amyloid 1-16	18.2	1956.1	0.1	
10	$\beta$ -amyloid 1-11	13.5	1326.3	0.001	
11	$\beta$ -conglycinin 165-178	5.2	1848.8	0.0001	

### 3-2. Digestion analyses




Fig 4. 1 pmol phosphorylase b Lys-C digestion fragments using ATHAP or CHCA.

Table 2. The Ion detection for 1 pmol phosphorylase b Lys-C digestion fragments and the SSRCalc Hydrophobicity<sup>a</sup>

No	phosphorylase b Lys-C digestion fragments		detection (+/-)	
	SSRCalc Hydrophobicity	m/z (Ave.)	ATHAP	CHCA
1	55.9	3602.2	++	-
2	53.9	3823.5	++	-
3	53.1	3890.3	++	-
4	51.0	3823.5	++	-
5	50.7	2198.6	++	+
6	45.8	2155.6	++	+
7	45.1	2742.0	+	+
8	42.8	2969.5	+	-
9	38.9	3504.9	++	-
10	35.1	1855.1	++	++
11	33.9	1657.0	+	+
12	33.6	2130.5	++	-
13	33.4	2629.0	++	++
14	31.4	1610.9	++	++
15	31.4	1814.1	++	++
16	31.2	1526.8	++	++
17	30.9	2043.3	++	++
18	30.4	1942.3	++	++
19	29.9	2449.7	-	+
20	28.4	1304.7	++	++
21	25.8	1263.4	-	++
22	24.7	1178.3	-	++
23	24.5	1290.5	+	++
24	18.9	1254.5	+	++
25	9.3	1102.2	-	++

- The **intensity was increased for hydrophobic peptides** with higher SSRCalc Hydrophobicity and **decreased for hydrophilic peptides** with lower SSRCalc Hydrophobicity **using ATHAP**.
- Hydrophobic peptide ions, which could not be detected by CHCA were detected by ATHAP**.
- The sequence coverage was increased by using the both results of ATHAP and CHCA.
  - 30% with CHCA or 49% with ATHAP
  - 55% with both the results of ATHAP and CHCA

<sup>a</sup> "++" indicates that the ions were detected with S/N  $\geq$  5, "+" indicates that the ions were detected with S/N = 2-5, and "-" indicates that the ions were not detected.

Table 3. Comparison of the Ion detection for 1 pmol phosphorylase b Lys-C digestion fragments using C6-, C8-, C10- or C12-ATHAP

no.	SSRCalc Hydrophobicity	phosphorylase b Lys-C digests				
		CHCA	C6-ATHAP	C8-ATHAP	C10-ATHAP	C12-ATHAP
1	55.9	-	-	++	++	-
2	53.9	-	-	++	++	-
3	53.1	-	-	++	+	-
4	51.0	-	-	++	++	-
5	50.7	+	+	++	++	-
6	45.8	+	+	++	++	-
7	45.1	+	-	+	+	-
8	42.8	-	-	+	+	-
9	38.9	-	-	++	++	-

- Hydrophobic peptide was detected with higher sensitivity using C8- (or C10-)ATHAP compared with C6-, C12-ATHAP or CHCA.

## 4. Conclusions

- A novel matrix **ATHAP** increased sensitivity of hydrophobic peptides by 10-fold.
- All the issues of ADHB were solved by ATHAP as follows:
  - The peptide ions were **detected throughout the entire matrix/analyte dried spot**.
  - ATHAP **works as a matrix itself**.
  - The detection of **hydrophilic peptide ions was suppressed**.

## References

- Fukuyama, Y.; Tanimura, R.; Maeda, K.; Watanabe, M.; Kawabata, S.; Iwamoto, S.; Izumi, S.; Tanaka, K. *Anal. Chem.* **2012**, *84*, 4237-4243.
- Fukuyama, Y.; Nakajima, C.; Furuichi, K.; Taniguchi, K.; Kawabata, S.; Izumi, S.; Tanaka, K. *Anal. Chem.* submitted.
- Krokhin, O.V.; *Anal. Chem.* **2006**, *78*, 7785-7795.

## Acknowledgments

This research is granted by the Japan Society for the Promotion of Science (JSPS) through the "Funding Program for World-Leading Innovative R&D on Science and Technology (FIRST Program)," initiated by the Council for Science and Technology Policy (CSTP).