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1. Abstract

Identification of bioactive peptides by DB searching is challenging because non-specific processing and PTMs results in very large search space. We developed novel sequence tag searching software to identify peptide variants effectively. For this purpose, focused tag library was built to store sequence tags obtained with high confidence from already identified peptides. As a preliminary test using small set of peptide data, correct sequence tags were derived effectively from undetermined spectra matching with tag library. We also confirmed correct peptides were identified just using top ranked sequence tag.

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Fig. 1. Functional/structural variety of endogenous peptide:
For example, pro-opiomelanocortin was cleaved enzymatically into various active peptides.

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
Fig. 2. Novel sequence tag search using focused peptide library

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
Fig. 3. Flow chart for process from tag query to tag database search:
 P_i : observed peptide ion mass, P_t : theoretical peptide ion mass matched with tag library.


Fig. 4. Schematic diagram of searching for PTM/substitution

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Table 1. Tryptic digest of actin for test 1

No.	Protein	Peptide	Precursor Ion Mass	Mod/Subst
1	β -actin (Human)	QEYDESGPSIVHR	1517 Da.	None
2	β -actin (Human)	QEYDESGPSIVHR	1545 Da.	N-term \rightarrow Formyl [27.8]
3	α -actin (Human)	QEYDEAGPSIVHR	1501 Da.	S \rightarrow A at Position 6

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