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Original Rev. 1.0 28/April/2013

Final Rev. 2.8 26/Feb/2015

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

Michael Karas, Franz Hillenkamp "Laser desorption ionization of proteins with molecular masses exceeding 10,000 daltons" Anal. Chem., Vol.60, October 1, p2299 (1988)

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

"SA" Sinapinic Acid 3,5-Dimethoxy-4-hydroxycinnamic acid	C₁₁H₁₂O₅ 530-59-6	224.06847 (224.210)		Ronald Beavis, Brian Chait "Cinnamic acid derivatives as matrices for ultraviolet laser desorption mass spectrometry of proteins" Rapid Commun. Mass Spectrom., Vol. 3, p432 (1989) Ronald Beavis, Brian Chait "Matrix-assisted laser-desorption mass spectrometry using 355 nm radiation" Rapid Commun. Mass Spectrom., Vol. 3, p436 (1989)
---	--	-------------------------------	--	---

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


Yutaka Aoki, Atsuhiko Toyama, Takashi Shimada, Tetsuyoshi Sugita, Chikage Aoki, Yukari Umino, Atsushi Suzuki, Daisuke Aoki, Yataro Daigo, Yusuke Nakamura, Taka-Aki Sato "A novel method for analyzing formalin-fixed paraffin embedded (FFPE) tissue sections by mass spectrometry imaging" *Proc. Jpn Acad. Ser.B*, Vol. 83, p205 (2007)

Hong Wang, Chee-Hong Wong, Alice Chin, Ayumu Taguchi, Allen Taylor, Samir Hanash, Sadanori Sekiya, Hidenori Takahashi, Masaki Murase, Shigeki Kajihara, Shinichi Iwamoto, Koichi Tanaka "Integrated mass spectrometry-based analysis of plasma glycoproteins and their glycan modifications" *Nature Protocols*, Vol. 6, p253 (2011)

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



Sadanori Sekiya, Yoshiki Yamaguchi, Koichi Kato, Koichi Tanaka "Mechanistic elucidation of the formation of reduced 2-aminopyridine-derivatized oligosaccharides and their application in matrix-assisted laser desorption/ionization mass spectrometry" *Rapid Commun. Mass Spectrom.*, Vol. 19, p3607 (2005)

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Yusaku Hioki, Ritsuko Tanimura, Shinichi Iwamoto, Koichi Tanaka "Nano-LC/MALDI-MS Using a Column-Integrated Spotting Probe for Analysis of Complex Biomolecule Samples" *Anal. Chem.*, Vol. 86, p2549 (2014)

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"CHCA" α-Cyano-4-hydroxy-cinnamic Acid	C₁₀H₇NO₃ 28166-41-8	189.04259 (189.167)		Ronald Beavis, T. Chaudhary, Brian Chait "α-Cyano-4-hydroxycinnamic acid as a matrix for matrix-assisted laser desorption mass spectrometry" Org. Mass Spectrom., Vol. 27, p156 (1992) Ole Vorm, Peter Roepstorff, Matthias Mann "Improved Resolution and Very High Sensitivity in MALDI TOF of Matrix Surfaces Made by Fast Evaporation" Anal. Chem., Vol. 66, p3281 (1994a)
--	--	-------------------------------	--	--

D.J.C. Pappin, P. Hojrup, A.J. Bleasby "Rapid identification of proteins by peptide-mass fingerprinting" Current Biol., Vol. 3, p327 (1993)


Toyofumi Nakanishi, Nobuhiko Okamoto, Koichi Tanaka, Akira Shimizu "Laser desorption time-of-flight mass spectrometric analysis of transferrin precipitated with antiserum: A unique simple method to identify molecular weight variants" Biol. Mass Spectrom., Vol. 23, p230 (1994b)

Plamen Demirev, Yen-Peng Ho, Victor Ryzhov, Catherine Fenselau "Microorganism Identification by Mass Spectrometry and Protein Database Searches" Anal. Chem., Vol. 71, p2732 (1999)


Yusaku Hioki, Hiroki Kuyama, Chikako Hamana, Kohei Takeyama, Koichi Tanaka "An improved sample preparation method for the sensitive detection of peptides by MALDI-MS" J. Mass Spectrom., Vol.48, p1217 (2013)

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

"MSA" 2-Methoxy-salicylic Acid 2-Hydroxy-5-methoxybenzoic Acid	C₈H₈O₄ 2612-02-4	168.04226 (168.147)		Todd Billeci, John Stults "Tryptic mapping of recombinant proteins by matrix-assisted laser desorption/ionization mass spectrometry" Anal. Chem., Vol. 65, p1709 (1993) Anthony Tsarbopoulos, Michael Karas, Kerstin Strupat, Birendra Pramanik, Tattanahalli Nagabhushan, Franz Hillenkamp "Comparative Mapping of Recombinant Proteins and Glycoproteins by Plasma Desorption and Matrix-Assisted Laser Desorption/Ionization Mass Spectrometry" Anal. Chem., Vol. 66, p2062 (1994) Arkady Gusev, William Wilkinson, Andrew Proctor, David Hercules "Improvement of signal reproducibility and matrix/comatrix effects in MALDI analysis" Anal. Chem., Vol. 67, p1034 (1995)
---	---	-------------------------------	--	--

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


Yuko Fukuyama, Shinichi Iwamoto, Koichi Tanaka "Rapid sequencing and disulfide mapping of peptides containing disulfide bonds by using 1,5-diaminonaphthalene as a reductive matrix" J. Mass Spectrom., Vol. 41, p191 (2006)


Yuko Fukuyama, Koichi Tanaka "Post-translational modification analyses using 1,5-diaminonaphthalene as a matrix" WP469 (ASMS2006)

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

Mass Spectra of Ribonuclease A by In Source Decay using 1,5-DAN (ASMS2006)

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Kratos/Shimadzu Axima TOF2

Uwe Pielers, Werner Zürcher, Martin Schär, Heinz Moser "Matrix-assisted laser desorption/ionization time-of-flight mass spectrometry: a powerful tool for the mass and sequence analysis of natural and modified oligonucleotides" Nucl. Acids Res., Vol. 21, p3191 (1993)

Nien-Yeen Hsu, Wen-Bin Yang, Chi-Huey Wong, Yuan-Chuan. Lee, Reiko T. Lee, Yi-Sheng Wang, Chung-Hsuan Chen "Matrix-assisted laser desorption/ionization mass spectrometry of polysaccharides with 2',4',6'-trihydroxyacetophenone as matrix" Rapid Commun. Mass Spectrom., Vol. 21, p2137 (2007a)

Gerald Stübiger, Omar Belgacem "Analysis of Lipids Using 2,4,6-Trihydroxyacetophenone as a Matrix for MALDI Mass Spectrometry" Anal. Chem., Vol. 79, p3206 (2007b)

PSD spectrum of $[M + Li]^+$ ions of sn-1-stearoyl-2-docosahexaenoyl-PC (22:6/18:0)

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Lucinda Cohen, Kerstin Strupat, Franz Hillenkamp "Analysis of quaternary protein ensembles by matrix assisted laser desorption/ionization mass spectrometry" J. Am. Soc. Mass Spectrom., Vol. 8, p1046 (1997)

Benoit Colsch, Amina Woods "Localization and imaging of sialylated glycosphingolipids in brain tissue sections by MALDI mass spectrometry" Glycobiol., Vol.20, p661 (2010)

"2,4-DHAP" 2',4'-Dihydroxyacetophenone	C₈H₈O₃ 89-84-9	152.04734 (152.147)		
--	---	-------------------------------	--	--

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


Fullerene

Salicylamide 2-hydroxybenzamide	C₇H₇NO₂ 65-45-2	137.04768 (137.136)		Kuang Jen Wu, Thomas Shaler, Christopher Becker "Time-of-flight mass spectrometry of underivatized single-stranded DNA oligomers by matrix-assisted laser desorption" Anal. Chem., Vol. 66, p1637 (1994)
---	--	-------------------------------	--	--

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Joern Krause, Markus Stoeckli, Urs Schlunegger "Studies on the selection of new matrices for ultraviolet matrix-assisted laser desorption/ionization time-of-flight mass spectrometry" Rapid Commun Mass Spectrom., Vol. 10, p1927 (1996)

Isovanillin 3-Hydroxy-4-methoxybenzaldehyde	C₈H₈O₃ 621-59-0	152.04734 (152.147)		Graeme Currie, John Yates "Analysis of oligodeoxynucleotides by negative-ion matrix-assisted laser desorption mass spectrometry" J. Am. Soc. Mass Spectrom., Vol. 4, p955 (1993)
---	--	-------------------------------	--	--

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Retinoic Acid all-trans-retinoic acid	C₂₀H₂₈O₂ 302-79-4	300.20893 (300.435)		Kathleen Walker, Michael Kahr, Charles Wilkins, Zhifu Xu, Jeffrey Moore "Analysis of hydrocarbon dendrimers by laser desorption time-of-flight and Fourier transform mass spectrometry" J. Am. Soc. Mass Spectrom., Vol. 5, p731 (1994)
---	--	-------------------------------	--	---

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David Schriemer, Liang Li "Detection of High Molecular Weight Narrow Polydisperse Polymers up to 1.5 Million Daltons by MALDI Mass Spectrometry" Anal. Chem., Vol. 68, p2721 (1996)

Anthranilic acid	C₇H₇NO₂ 118-92-3	137.04768 (137.136)		Bing Wang, Klaus Biemann "Matrix-Assisted Laser Desorption / Ionization Time-of-Flight Mass Spectrometry of Chemically Modified Oligonucleotides" Anal. Chem., Vol. 66, p1918 (1994)
-------------------------	---	-------------------------------	--	--

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
“3-AQ/CHCA” (3-AQ, CHCA説明を参照)


Kumar Kolli, Ron Orlando “A new matrix for MALDI on magnetic sector instrument with point detectors” Rapid Commun. Mass Spectrom., Vol. 10, p923 (1996)

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
Kaoru Kaneshiro, Yuko Fukuyama, Shinichi Iwamoto, Sadanori Sekiya, Koichi Tanaka “Highly Sensitive MALDI Analyses of Glycans by a New Aminoquinoline-Labeling Method Using 3-Aminoquinoline / α -Cyano-4-hydroxycinnamic Acid Liquid Matrix” Anal. Chem., Vol. 83, p3663 (2011)
Kaoru Kaneshiro, Makoto Watanabe, Kazuya Terasawa, Hiromasa Uchimura, Yuko Fukuyama, Shinichi Iwamoto, Taka-Aki Sato, Kazuharu Shimizu, Koichi Tanaka “Rapid Quantitative Profiling of N-Glycan by the Glycan-Labeling Method Using 3-Aminoquinoline/ α -Cyano-4-hydroxycinnamic Acid” Anal. Chem., Vol. 84, p7146 (2012a)

Koichi Tanaka, Sadanori Sekiya, Masafumi Jinno, Makoto Hazama, Kei Kodera, Shinichi Iwamoto “Macromolecule Measurement using MALDI-DIT-MS” MP011 (ASMS2008)


Koichi Tanaka, Sadanori Sekiya, Shinichi Iwamoto “MALDI-DITMS/MS for High Mass High Sensitivity and High Resolution Measurement” MP586 (ASMS2009)

Sadanori Sekiya, Ken-ichi Taniguchi, Koichi Tanaka “On-target separation of analyte with 3-aminoquinoline / α -cyano-4-hydroxycinnamic acid liquid matrix for matrix-assisted laser desorption / ionization mass spectrometry” Rapid Commun. Mass Spectrom., Vol. 26, p693 (2012c)

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

Takuya Shuo, Naohiko Koshikawa, Daisuke Hoshino, Tomoko Minegishi, Hiroko Ao-Kondo, Masaaki Oyama, Sadanori Sekiya, Shinichi Iwamoto, Koichi Tanaka, Motoharu Seiki “Detection of the Heterogeneous O-Glycosylation Profile of MT1-MMP Expressed in Cancer Cells by a Simple MALDI-MS Method” PLOS ONE, Vol. 7, e43751 (2012d)

Takashi Nishikaze, Kaoru Kaneshiro, Shin-ichiro Kawabata, Koichi Tanaka “Structural Analysis of N-Glycans by the Glycan-Labeling Method Using 3-Aminoquinoline-Based Liquid Matrix in Negative-Ion MALDIMS” Anal. Chem., Vol. 84, p9453 (2012b)

Acquired in 2009 using Axima-TOF² (Unpublished)

Yuko Fukuyama, Kohei Takeyama, Shin-ichiro Kawabata, Shinichi Iwamoto, Koichi Tanaka “An optimized matrix-assisted laser desorption/ionization sample preparation using a liquid matrix, 3-aminoquinoline/ α -cyano-4-hydroxycinnamic acid, for phosphopeptides” Rapid Commun. Mass Spectrom., Vol. 26, p2454 (2012e)

Makoto Watanabe, Kazuya Terasawa, Kaoru Kaneshiro, Hiromasa Uchimura, Rie Yamamoto, Yuko Fukuyama, Kazuharu Shimizu, Taka-Aki Sato, Koichi Tanaka “Improvement of mass spectrometry analysis of glycoproteins by MALDI-MS using 3-aminoquinoline/ α -cyano-4-hydroxycinnamic acid” Anal. Bioanal. Chem., Vol. 405, p4289 (2013)

Takashi Nishikaze, Shin-ichiro Kawabata, Koichi Tanaka “Fragmentation Characteristics of Deprotonated N-linked Glycopeptides: Influences of Amino Acid Composition and Sequence” J. Am. Soc. Mass Spectrom., Vol. 25, p988 (2014)

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Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹⁾	Structure	Reference(s)
“CSA” 5-Chloro-salicylic Acid 5-Chloro-2-hydroxybenzoic acid	C₇H₅ClO₃ 321-14-2	171.993 (172.566)		Lisa Preston-Schaffter, Gary Kinsel, David Russell “Effects of heavy-atom substituents on matrices used for matrix-assisted laser desorption/ionization mass spectrometry” J. Am. Soc. Mass Spectrom., Vol. 5, p800 (1994)
DHB 5位 OHがClに置換した構造。疎水性化合物(非極性Polymer、膜タンパク質)測定等が報告されている。				
2-bromo-4,6-dinitroaniline	C₆H₄BrN₃O₄ 1817-73-8	260.93852 (262.018)		Burkhard Rosinke, Kerstin Strupat, Franz Hillenkamp, Jürg Rosenbusch, Norbert Dencher, Ulrike Krüger, Hans-Joachim Galla “Matrix-assisted laser desorption/ionization mass spectrometry (MALDI-MS) of membrane proteins and non-covalent complexes” J. Mass Spectrom., Vol. 30, p1462 (1995)
5-bromo-2,4-dihydroxybenzoic acid (monohydrate)	C₇H₅BrO₄ 7355-22-8	231.93712 (233.016)		Josef Chmelík, Josef Planeta, Pavel Řehulka, Josef Chmelík “Determination of molecular mass distribution of silicone oils by supercritical fluid chromatography, matrix-assisted laser desorption/ionization time-of-flight mass spectrometry and their off-line combination” J. Mass Spectrom., Vol. 36, p760 (2001)
4-Hydroxy-3-iodo-5-nitrobenzoic acid	C₇H₄INO₅ 10463-17-9	308.91342 (309.014)		D. J. Harvey “Matrix-assisted laser desorption/ionization mass spectrometry of phospholipids” J. Mass Spectrom., Vol. 30, p1333 (1995)
従来Matrix構造(の水素部位)にCl, Br, I等 heavy-atomを導入、excited state relaxationを変化させ、イオン化効率を高める事等が行われている。				
1,4-Dihydroxy-2-naphthoic acid	C₁₁H₈O₄ 31519-22-9	204.04226 (204.179)		Xiaodong Tang, Peter Dreifuss, Akos Vertes “New matrices and accelerating voltage effects in matrix-assisted laser desorption/ionization of synthetic polymers” Rapid Commun. Mass Spectrom., Vol. 9, p1141 (1995)
合成高分子・脂質等の測定に試行された。				
“HIQ” 1-hydroxy isoquinoline	C₉H₇NO 491-30-5	145.05276 (145.158)		Martin Mohr, K. OlafBörnsen, H. Michael Widmer “Matrix-assisted laser desorption/ionization mass spectrometry: Improved matrix for oligosaccharides” Rapid Commun. Mass Spectrom., Vol. 9, p809 (1995)
(DHB等に1~数10%添加する事により)bufferやその他contaminantが存在しても、糖鎖やタンパク質をquality高く測定できる。				
Indole-2-carboxylic acid	C₉H₇NO₂ 1477-50-5	161.04768 (161.157)		Hanno Ehring, Bo Sundqvist “Studies of the MALDI process by luminescence spectroscopy” J. Mass Spectrom., Vol. 30, p1303 (1995)
IR-Laserでも用いられる。				
Graphite for “Graphite Assisted Laser Desorption / Ionization” (GALDI)				Jan. Sunner, Edward. Dratz, Yu-Chie. Chen “Graphite surface-assisted laser desorption/ ionization time-of-flight mass spectrometry of peptides and proteins from liquid solutions” Anal. Chem., Vol. 67, p4335 (1995)
Michael Dale, Richard Knochenmuss, Renato Zenobi “Graphite / Liquid Mixed Matrices for Laser Desorption / Ionization Mass Spectrometry” Anal. Chem., Vol. 68, p3321 (1996)				
Stefan Zumbuhl, Richard Knochenmuss Stefan Wulfert, Frederic Dubois, Michael Dale, Renato Zenobi “A Graphite-Assisted Laser Desorption/Ionization Study of Light-Induced Aging in Triterpene Damar and Mastic Varnishes” Anal. Chem., Vol. 70, p707 (1998)				
Claudie Black, Chris Poile, John Langley, Julie Herniman “The use of pencil lead as a matrix and calibrant for matrix-assisted laser desorption/ionisation” Rapid Commun. Mass Spectrom., Vol. 20, p1053 (2006)				

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Polymer-Assisted Laser Desorption / Ionization (PALDI)用Matrixの一種。

"Proton sponge" "DMAN" N,N,N',N'-Tetramethyl-1,8-naphthalenediamine	C₁₄H₁₈N₂ 20734-58-1	214.14700 (214.306)		James Stephenson Jr., Scott McLuckey "Ion/ion reactions in the gas phase: Proton transfer reactions involving multiply-charged proteins" J. Am. Chem. Soc., Vol. 118, p7390 (1996) Rohit Shroff, Ales Svatos "Proton Sponge: A Novel and Versatile MALDI Matrix for the Analysis of Metabolites Using Mass Spectrometry" Anal. Chem., Vol.81, p7954 (2009)
---	--	-------------------------------	--	---

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Christina Siedschlag, Heinrich Luftmann, Christian Wolff, Jochen Mattay "Functionalization of [60]fullerene by photoinduced electron transfer (PET): Syntheses of 1-substituted 1,2-dihydro[60]fullerenes" Tetrahedron, Vol. 53, p3587 (1997)



Mark Wyatt, Bridget Stein, Gareth Brenton "Characterization of Various Analytes Using Matrix-Assisted Laser Desorption/Ionization Time-of-Flight Mass Spectrometry and 2-[(2E)-3-(4-tert-Butylphenyl)-2-methylprop-2-enylidene] malononitrile Matrix" Anal. Chem., Vol. 78, p199 (2006)

Shin-ichirou Kawabata, Shinichi Iwamoto, Koichi Tanaka "Stable Ion Formation Techniques with MALDI and Their Applications on High Mass Substances" 16P-039 (MSSJ2010)

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

5-nitrosalicylic acid	C₇H₅NO₅ 96-97-9	183.01677 (183.118)		J. Guittard, J. C. Blais, G. Bolbach, A. Brunot, M. Tessier, E. Marechal, J. C. Tabet "Ionization of synthetic polymers in matrix-assisted laser desorption ionization (MALDI)" Proc. 45th ASMS Conf. Mass Spectrom. Allied Topics, p840 (1997) Daiki Asakawa, Mitsuo Takayama "Specific cleavage at peptide backbone Cα-C and CO-N bonds during matrix-assisted laser desorption/ionization in-source decay mass spectrometry with 5-nitrosalicylic acid as the matrix" Rapid Commun. Mass Spectrom., Vol. 25, p2379 (2011)
------------------------------	--	-------------------------------	--	---

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Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹	Structure	Reference(s)
nor-harman 9H-pyrido[3,4-b]indole	C₁₁H₈N₂ 244-63-3	168.06875 (168.195)		Hiroshi Nonami, Shinsaku Fukui, Rosa Erra-Balsells "β-Carboline Alkaloids as Matrices for Matrix-assisted Ultraviolet Laser Desorption Time-of-flight Mass Spectrometry of Proteins and Sulfated Oligosaccharides: a Comparative Study Using Phenylcarbonyl Compounds, Carbazoles and Classical Matrices" J. Mass Spectrom., Vol.32, p287 (1997)
Harmane	C₁₂H₁₀N₂ 486-84-0	182.08440 (182.221)		Koichi Tanaka, Hiroshi Nonami, Yuko Fukuyama, Rosa Erra-Balsells "β-Carbolines as Matrices for MALDI ToF-MS(MS) in Positive and Negative Modes" WP222 (ASMS1998)
Harmine	C₁₃H₁₂N₂O 442-51-3	212.09496 (212.247)		Hiroshi Nonami, Koichi Tanaka, Yuko Fukuyama, Rosa Erra-Balsells "β-Carboline Alkaloids as Matrices for UV-Matrix assisted Laser Desorption/Ionization Time-of-flight Mass Spectrometry in Positive and Negative Ion Modes. Analysis of Proteins of High Molecular Mass, and of Cyclic and Acyclic Oligosaccharides" Rapid Commun. Mass Spectrom., Vol. 12, p285 (1998)
Harmol	C₁₂H₁₀N₂O 487-03-6	198.07931 (198.221)		Hiroshi Nonami, Feiyue Wu, Randolph Thummel, Yuko Fukuyama, Hidenobu Yamaoka, Rosa Erra-Balsells "Evaluation of pyridoindoles, pyridylindoles and pyridylpyridoindoles as matrices for ultraviolet matrix-assisted laser desorption/ionization time-of-flight mass spectrometry" Rapid Commun. Mass Spectrom. Vol. 15, p2354 (2001)
Harmaline	C₁₃H₁₄N₂O 304-21-2	214.11061 (214.263)		Mariana Barboza, Vilma Duschak, Yuko Fukuyama, Hiroshi Nonami, Rosa Erra-Balsells, Juan Cazzulo, Alicia Couto "Structural analysis of the N-glycans of the major cysteine proteinase of Trypanosoma cruzi -- Identification of sulfated high-mannose type oligosaccharides --" FEBS Journal, Vol. 272, p3803 (2005)
Harmalol	C₁₂H₁₂N₂O 525-57-5	200.09496 (200.236)		"MALDI-MS Technical Reports" No.08

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

Ref: ΔpKa = [pKa(Lower Singlet Excited State)] - [pKa(Ground State)]

nor-Harmane Solubilities / Stabilities (>1 Month) to Several Solvents

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
Only these compounds that possessed indole N-H and pyridine nitrogen groups in a 1-4 relationship demonstrated good matrix properties with a nitrogen laser. (2001論文)

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


“CHCAB”			$\left[\begin{array}{c} \text{H} \\ \\ \text{C}_4\text{H}_9\text{-N}^+\text{-H} \\ \\ \text{H} \end{array} \right] \left[\begin{array}{c} \text{CH=C(CN)COO}^- \\ \\ \text{OH} \end{array} \right]$	(上記論文に加え 下記も) Yuko Fukuyama, Natsumi Funakoshi, Shinichi Iwamoto, Koichi Tanaka “Adding methanol to α-cyano-4-hydroxycinnamic acid butylamine salt as a liquid matrix to form a homogeneous spot on a focusing plate for highly sensitive and reproducible analyses in matrix-assisted laser desorption/ionization mass spectrometry” <i>Rapid Commun. Mass Spectrom.</i> , Vol. 28, p662 (2014)
----------------	--	--	---	---

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Takahiro Harada, Akiko Yuba-Kubo, Yuki Sugiura, Nobuhiro Zaima, Takahiro Hayasaka, Naoko Goto-Inoue, Masatoshi Wakui, Makoto Suematsu, Kengo Takeshita, Kiyoshi Ogawa, Yoshikazu Yoshida, Mitsutoshi Setou "Visualization of Volatile Substances in Different Organelles with an Atmospheric-Pressure Mass Microscope" Anal. Chem., Vol. 81, p9153 (2009)


Akiko Kubo, Mitsuyo Ohmura, Masatoshi Wakui, Takahiro Harada, Shigeki Kajihara, Kiyoshi Ogawa, Hiroshi Suemizu, Masato Nakamura, Mitsutoshi Setou, Makoto Suematsu "Semi-quantitative analyses of metabolic systems of human colon cancer metastatic xenografts in livers of superimmunodeficient NOG mice" Anal. Bioanal. Chem., Vol. 400, p1895 (2011)

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
"iMScope TRIO"


Imaging MALDI (Prototype)
(下記データはPrototype機で測定)

Masahiro Kawashima, Noriko Iwamoto, Nobuko Kawaguchi-Sakita, Masahiro Sugimoto, Takayuki Ueno, Yoshiaki Mikami, Kazuya Terasawa, Taka-Aki Sato, Koichi Tanaka, Kazuharu Shimizu, Masakazu Toi "High-Resolution Imaging Mass Spectrometry Reveals Detailed Spatial Distribution of Phosphatidylinositols in Human Breast Cancer" Cancer Sci., Vol. 104, p1372 (2013)

Takayuki Goto, Naoki Terada, Takahiro Inoue, Kenji Nakayama, Yoshiyuki Okada, Takeshi Yoshikawa, Yu Miyazaki, Masayuki Uegaki, Shinji Sumiyoshi, Takashi Kobayashi, Tomomi Kamba, Koji Yoshimura, Osamu Ogawa "The Expression Profile of Phosphatidylinositol in High Spatial Resolution Imaging Mass Spectrometry as A Potential Biomarker for Prostate Cancer" PLOS ONE, Vol. 9, e90242 (2014)

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Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹	Structure	Reference(s)
“DNPH” 2,4-Dinitrophenyl hydrazine	C₆H₆N₄O₄ 119-26-6	198.03890 (198.136)		François Fenaïlle, Jean-Claude Tabet, Philippe Guy “Identification of 4-Hydroxy-2-nonenal-Modified Peptides within Unfractionated Digests Using Matrix-Assisted Laser Desorption/Ionization Time-of-Flight Mass Spectrometry” <i>Anal. Chem.</i> , Vol. 76, p867 (2004) Kristin Teuber, Maria Fedorova, Ralf Hoffmann, Jürgen Schiller “2,4-Dinitrophenylhydrazine as a New Reactive Matrix to Analyze Oxidized Phospholipids by MALDI-TOF Mass Spectrometry” <i>Analytical Letters</i> Vol. 45, p968 (2012)

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
“DHB/B” 2,5-Dihydroxybenzoic Acid Butylamine				Marko Mank, Bernd Stahl, Günther Boehm “2,5-Dihydroxybenzoic Acid Butylamine and Other Ionic Liquid Matrixes for Enhanced MALDI-MS Analysis of Biomolecules” <i>Anal. Chem.</i> , Vol. 76, p2938 (2004) Masoud Zabet-Moghaddam, Elmar Heinze, Andreas Tholey “Qualitative and quantitative analysis of low molecular weight compounds by ultraviolet matrix-assisted laser desorption/ionization mass spectrometry using ionic liquid matrices” <i>Rapid Commun Mass Spectrom</i> , Vol. 18, p141 (2004)
--	--	--	--	---

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Pyrene	C₁₆H₁₀ 129-00-0	202.07825 (202.251)		Stephen Macha, Patrick Limbach, Philip Savickas “Application of Nonpolar Matrices for the Analysis of Low Molecular Weight Nonpolar Synthetic Polymers by Matrix-Assisted Laser Desorption/Ionization Time-of-Flight Mass Spectrometry” <i>J. Am. Soc. Mass Spectrom.</i> , Vol. 11, p731 (2000) Junko Amano, Koichi Tanaka “Negative-ion MALDI MS ⁿ for Identification and Quantification of Isomeric Fucosylated Oligosaccharide Mixtures” <i>ThP026 (ASMS2005)</i> Junko Amano, Fumio Tougasakoi, Koichi Tanaka “Highly sensitive MALDI MS for Identification of Polysulfated Oligosaccharides by a Simple Pyrene-derivatization Method” <i>ThP298 (ASMS2006)</i> Junko Amano, Daisuke Sugahara, Kenji Osumi Koichi Tanaka “Negative-ion MALDI-QIT-TOFMS ⁿ for structural determination of fucosylated and sialylated oligosaccharides labeled with a pyrene derivative” <i>Glycobiol.</i> , Vol. 19, p592 (2009) Junko Amano, Minako Osanai, Takahiro Orita, Daisuke Sugahara, Kenji Osumi “Structural determination by negative-ion MALDI-QIT-TOFMS ⁿ after pyrene derivatization of variously fucosylated oligosaccharides with branched decaose cores from human milk” <i>Glycobiol.</i> Vol. 19, p601 (2009) Junko Amano, Takashi Nishikaze, Fumio Tougasaki, Hiroshi Jinmei, Ichiro Sugimoto, Shu-ichi Sugawara, Masaya Fujita, Kenji Osumi, Mamoru Mizuno “Derivatization with 1-Pyrenyldiazomethane Enhances Ionization of Glycopeptides but Not Peptides in Matrix-Assisted Laser Desorption/Ionization Mass Spectrometry” <i>Anal. Chem.</i> , Vol. 82, p8738 (2010)
---------------	---	-------------------------------	--	---

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
Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹	Structure	Reference(s)
“DEA/CHCA” α-Cyano-4-hydroxycinnamic acid diethylammonium salt diethylammonium α-cyano-4-hydroxy cinnamate				Ying Li, Michael Gross “Ionic-liquid matrices for quantitative analysis by MALDI-TOF mass spectrometry” J. Am. Soc. Mass Spectrom., Vol. 15, p1833 (2004)

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“3H4NBA” 3-hydroxy-4-nitrobenzoic acid	C₇H₅NO₅ 619-14-7	183.01677 (183.118)	
--	---	-------------------------------	--

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Ei-ichi Matsuo, Chikako Toda, Makoto Watanabe, Noriyuki Ojima, Shunsuke Izumi, Koichi Tanaka, Susumu Tsunasawa, Osamu Nishimura “Selective detection of 2-nitrobenzenesulfonyl-labeled peptides by matrix-assisted laser desorption / ionization time of flight mass spectrometry using a novel matrix” Proteomics Vol. 6, p2042 (2006)

Watanabe M, Takemasa I, Nishimura N, Matsubara T, Yoshioka S, Miyake M, Nagai K, Monden M, Nishimura O “The proteomic profiling of human colorectal cancer tissue using NBS method” MP624 (ASMS2006)

Keli Ou, Djohan Kesuma, Kumaresan Ganesan, Kun Yu, Sou Yen Soon, Suet Ying Lee, Xin Pei Goh, Michelle Hooi, Wei Chen, Hiroyuki Jikuya, Tetsuo Ichikawa, Hiroki Kuyama, Ei-ichi Matsuo, Osamu Nishimura, Patrick Tan “Quantitative Profiling of Drug-Associated Proteomic Alterations by Combined 2-Nitrobenzenesulfonyl Chloride (NBS) Isotope Labeling and 2DE/MS Identification” J. Proteome Res., Vol. 5, p2194 (2006)


Ueda K, Katagiri T, Shimada T, Irie S, Sato T, Nakamura Y, Daigo Y “Comparative profiling of serum glycoproteome by sequential purification of glycoproteins and 2-nitrobenzenesulfonyl (NBS) stable isotope labeling: a new approach for the novel biomarker discovery for cancer” J Proteome Res., Vol. 6, p3475 (2007)

Okamura N, Masuda T, Gotoh A, Shirakawa T, Terao S, Kaneko N, Suganuma K, Watanabe M, Matsubara T, Seto R, Matsumoto J, Kawakami M, Yamamori M, Nakamura T, Yagami T, Sakaeda T, Fujisawa M, Nishimura O, Okumura K. “Quantitative proteomic analysis to discover potential diagnostic markers and therapeutic targets in human renal cell carcinoma” Proteom., Vol. 8, p3194 (2008)

Watanabe M, Takemasa I, Kawaguchi N, Miyake M, Nishimura N, Matsubara T, Matsuo E, Sekimoto M, Nagai K, Matsuura N, Monden M, Nishimura O “An application of the 2-nitrobenzenesulfonyl method to proteomic profiling of human colorectal carcinoma: A novel approach for biomarker discovery” Proteomics Clin Appl., Vol. 2, p925 (2008)

Ei-ichi Matsuo, Makoto Watanabe, Hiroki Kuyama, Osamu Nishimura “A new strategy for protein biomarker discovery utilizing 2-nitrobenzenesulfonyl (NBS) reagent and its applications to clinical samples” J. Chromatogr. B, Vol. 877, p2607 (2009)

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Yu Fu, Songyun Xu, Chensong Pan, Mingliang Ye, Hanfa Zou, Baochuan Guo “A matrix of 3,4-diaminobenzophenone for the analysis of oligonucleotides by matrix-assisted laser desorption/ionization time-of-flight mass spectrometry” Nucl. Acids Res., Vol. 34, e94 (2006)

Songyun Xu, Mingliang Ye, Danke Xu, Xin Li, Chensong Pan, Hanfa Zou “Matrix with High Salt Tolerance for the Analysis of Peptide and Protein Samples by Desorption / Ionization Time-of-Flight Mass Spectrometry” Anal. Chem., Vol. 78, p2593 (2006)

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Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹	Structure	Reference(s)
“Nano-PALDI” Nanoparticle-Assisted Laser Desorption / Ionization				Shu Taira, Yuki Sugiura, Shinji Moritake, Shuichi Shimma, Yuko Ichiyanagi “Nanoparticle-assisted laser desorption/ionization for ultra resolution imaging mass spectrometry” WP068 (ASMS 2007) Shu Taira, Yuki Sugiura, Shinji Moritake, Shuichi Shimma, Yuko Ichiyanagi, Mitsutoshi Setou “Nanoparticle-Assisted Laser Desorption / Ionization Based Mass Imaging with Cellular Resolution” Anal. Chem., Vol. 80, p4761 (2008)

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
“CA” <i>trans-p-Coumaric Acid</i> <i>trans-4-hydroxycinnamic acid</i>	C₉H₈O₃ 501-98-4	164.04734 (164.158)		
--	--	-------------------------------	--	--

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

“G₂/CHCA”				Tatiana Laremore, Fuming Zhang, Robert Linhardt “Ionic Liquid Matrix for Direct UV-MALDI-TOF-MS Analysis of Dermatan Sulfate and Chondroitin Sulfate Oligosaccharides” Anal. Chem., Vol. 79, p1604 (2007)
“G₃/CA”				Yuko Fukuyama, Shuichi Nakaya, Yuzo Yamazaki, Koichi Tanaka “Improvement of sulfated / sialylated / neutral oligosaccharide and glycopeptide analyses using ionic liquid matrices” WP290 (ASMS 2007) Yuko Fukuyama, Atsuhiko Toyama, Koichi Tanaka “On-target digestion using liquid matrixes for glycopeptide analysis” ThP382 (ASMS 2008) Yuko Fukuyama, Shuichi Nakaya, Yuzo Yamazaki, Koichi Tanaka “Ionic Liquid Matrixes Optimized for MALDI-MS of Sulfated / Sialylated / Neutral Oligosaccharides and Glycopeptides” Anal. Chem., Vol. 80, p2171 (2008) Takashi Nishikaze, Yuko Fukuyama, Shin-ichirou Kawabata, Koichi Tanaka “Sensitive Analyses of Neutral N-Glycans using Anion-Doped Liquid Matrix G ₃ CA by Negative-Ion Matrix-Assisted Laser Desorption/Ionization Mass Spectrometry” Anal. Chem., Vol. 84, p6097 (2012) Takashi Nishikaze, Shin-ichirou Kawabata, Shinichi Iwamoto and Koichi Tanaka “Reversible hydrazide chemistry-based enrichment for O-GlcNAc-modified peptides and glycopeptides having non-reducing GlcNAc residues” Analyst, Vol. 138, p7224 (2013)

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
Negative-ion mass spectra of Ribonuclease B digests with DHB and G₃CA

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Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹⁾	Structure	Reference(s)
Aminopyrazine	C₄H₅N₃ 5049-61-6	95.04835 (95.103)		Muhammad Hashir, Guenther Stecher, Guenther Bonn "Identification of low molecular weight carbohydrates employing new binary mixtures for matrix-assisted laser desorption/ionisation mass spectrometry" Rapid Commun. Mass Spectrom., Vol. 22, p2185 (2008)

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"IMTBA/CHCA" N-isopropyl-N-methyl-N-tert-butylammonium α-cyano-4-hydroxycinnamate				Jeffrey Crank, Daniel Armstrong "Towards a Second Generation of Ionic Liquid Matrices (ILMs) for MALDI-MS of Peptides, Proteins, and Carbohydrates" J. Am. Soc. Mass Spectrom., Vol. 20, p1790 (2009) Carlos Serrano, Yi Zhang, Jian Yang, Kevin Schug "Matrix - assisted laser desorption / ionization mass spectrometric analysis of aliphatic biodegradable photoluminescent polymers using new ionic liquid matrices" Rapid Commun. Mass Spectrom., Vol. 25, p1152 (2011)
"DIEA/CHCA" N,N-diisopropylethylammonium α-cyano-4-hydroxycinnamate				
"DIEA/FA" N,N-diisopropylethylammonium ferulate				

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

"5-ASA" 5-amino-salicylic acid	C₇H₇NO₃ 89-57-6	153.04259 (153.135)		Motoshi Sakakura, Mitsuo Takayama "In-source decay and fragmentation characteristics of peptides using 5-aminosalicylic acid as a matrix in matrix-assisted laser desorption/ionization mass spectrometry" J. Am. Soc. Mass Spectrom., Vol. 21, p979 (2010)
--	--	-------------------------------	--	---

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MALDI-MS spectra of [pTyr6]-ACTH18-35 obtained with 1,5-DAN (top) and 5-ASA(bottom)

Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹⁾	Structure	Reference(s)
5-aryl-1H-tetrazoles	----- -----		<p>R: H, 2-Br, 4-OH</p>	<p>m: matrix M: Propranolol HCl</p> <p>CHCA</p> <p>DHB</p> <p>5-phenyl-1H-tetrazole</p>
<p>中越雅道, 高橋豊, 船越なつ美, 中島宏樹 “低分子化合物測定用新規MALDIマトリックスの開発(1)” 2P-038 (MSSJ2009)</p> <p>中島宏樹, 船越なつ美, 岩本慎一, 高橋豊, 中越雅道 “低分子化合物測定用新規MALDIマトリックスの開発(2)~低分子化合物の測定~” 2P-039 (MSSJ2009)</p> <p>Masamichi Nakakoshi, Natsumi Funakoshi, Hiroki Nakajima, Masayoshi Tsubuki “New Type Of Matrix For Analysis Of Low Molecular Weight Compounds by MALDI-TOF MS” HETEROCYCLES, Vol. 85, p927 (2012)</p> <p>Matrix関連ピークが検出され難いため、低分子化合物の測定に有用、と考えられる。</p>				


“DiFCCA” α-Cyano-2,4-difluorocinnamic acid	C₁₀H₅F₂NO₂ 953905-51-6	209.02884 (209.148)		Kristin Teuber, Jürgen Schiller, Beate Fuchs, Michael Karas, Thorsten Jaskolla “Significant sensitivity improvements by matrix optimization: a MALDI-TOF mass spectrometric study of lipids from hen egg yolk” Chemistry and Physics of Lipids, Vol. 163, p552 (2010) Jens Soltwisch, Thorsten Jaskolla, Klaus Dreisewerd “Color Matters—Material Ejection and Ion Yields in UV-MALDI Mass Spectrometry as a Function of Laser Wavelength and Laser Fluence” J Am. Soc. Mass Spectrom., Vol. 24, p1477 (2013)
LipidsをPositive modeで測定する場合、DiFCCAは他のMatrixよりも1桁以上感度が高いが、Negative modeでは特にphospholipids測定に不向きである。				

3-Nitrobenzonitrile	C₇H₄N₂O₂ 619-24-9	148.02728 (148.119)		Shirley Lomeli, Ivory Peng, Sheng Yin, Rachel Loo, Joseph Loo “New Reagents for Increasing ESI Multiple Charging of Proteins and Protein Complexes” J Am. Soc. Mass Spectrom., Vol. 21, p127 (2010) Ellen Inutan, Sarah Trimpin “Matrix assisted ionization vacuum (MAIV), a new ionization method for biological materials analysis using mass spectrometry” Mol. Cell. Proteom., Vol. 12, p792 (2010)
主に Matrix Assisted Ionization Vacuum (MAIV) 用Matrixとして用いられる。				


“ADHB” O-alkylated dihydroxy benzoic acid	C₁₅H₂₂O₄ -----	266.15181 (266.333)	<p>DHB structure</p> <p>Hydrophobic alkyl group</p> <p>構造修正論文: Anal. Chem., Vol. 86, p 5187 (2014)</p>	福山裕子, 谷村里都子, 泉俊輔, 岩本慎一, 田中耕一 “新規添加剤ADHB混合マトリックスを用いた疎水性ペプチド高感度分析” 質量分析総合討論会 2P-40 (MSSJ 2011) Yuko Fukuyama, Ritsuko Tanimura, Kazuki Maeda, Makoto Watanabe, Shin-ichirou Kawabata, Shinichi Iwamoto, Shunsuke Izumi, Koichi Tanaka “Alkylated Dihydroxybenzoic Acid as a MALDI Matrix Additive for Hydrophobic Peptide Analysis” Anal. Chem., Vol. 84, p4237 (2012)
---	---	-------------------------------	--	--

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Name	Empir. Form. CAS No.	Monoiso. Mass (Aver. Mass) ¹	Structure	Reference(s)
------	----------------------	---	-----------	--------------

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

Yuko Fukuyama, Natsumi Funakoshi, Kohei Takeyama, Yusaku Hioki, Takashi Nishikaze, Kaoru Kaneshiro, Shin-ichiro Kawabata, Shinichi Iwamoto, Koichi Tanaka "3-Aminoquinoline / *p*-Coumaric Acid as a MALDI Matrix for Glycopeptides, Carbohydrates, and Phosphopeptides" Anal. Chem., Vol. 86, p1937 (2014)

Takashi Nishikaze, Shin-ichiro Kawabata, Koichi Tanaka "In-Depth Structural Characterization of N-Linked Glycopeptides Using Complete Derivatization for Carboxyl Groups Followed by Positive- and Negative-Ion Tandem Mass Spectrometry" Anal. Chem., Vol. 86, p5360 (2014)

Ritsuko Yoda, Yusaku Hioki, Takashi Nishikaze, Shinichi Iwamoto, Koichi Tanaka "Sequential detection of glycopeptides and their glycans by using LC/MALDI-MS and on-target digestion" SCS25, P-06 (2014)

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

“ATHAP”

Alkylated Trihydroxyacetophenone



252.13616 (252.306)

Yuko Fukuyama, Chihiro Nakajima, Keiko Furuichi, Kenichi Taniguchi, Shin-ichiro Kawabata, Shunsuke Izumi, Koichi Tanaka "Alkylated Trihydroxyacetophenone as a MALDI Matrix for Hydrophobic Peptides" Anal. Chem., Vol. 85, p9444 (2013)

Yuko Fukuyama, Chihiro Nakajima, Shunsuke Izumi, Koichi Tanaka "Enhanced MALDI-MS using Alkylated Trihydroxyacetophenone as a Matrix for Hydrophobic Peptides in Membrane Protein Analysis" MP130 (ASMS2014)

Mass spectra of Intact Bacteriorhodopsin(BR)

Mass spectra of BR trypsin digestion (100 fmol)
★: Hydrophobic peptides containing TM domains
●: Hydrophilic peptides

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